# B4<sup>®</sup>

DENTSPLY

## Pre-Impression Surface Optimizer

### DIRECTIONS FOR USE – ENGLISH

For dental use only.  
USA: Rx only.

B4<sup>®</sup> Pre-Impression Surface Optimizer is designed to be used as a pre-impresion conditioning solution. When applied to substrates such as sulcular tissue, tooth preparations, ceramic, composite or metal core build-up material or implant abutments or copings, it provides a uniform surface condition with reduced surface tension prior to the application of Aquasil Ultra Impression Material. This standardized surface with reduced contact angle on all impressed substrates allows capture of fine preparation detail both supra and subgingivally.

### COMPOSITION

Polyethylene Glycol (PEG); Ethanol; Water; Surfactant; Pigments

### INDICATIONS

Preparation conditioning solution for uniform surface wetting prior to elastomeric impression procedures.

### CONTRAINDICATIONS

B4<sup>®</sup> Optimizer is contraindicated for direct application to dental pulp tissue.

#### Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

### WARNINGS

1. B4<sup>®</sup> Optimizer contains PEG (polyethylene glycol). Avoid prolonged or repeated contact with skin (allergic contact dermatitis), oral soft tissues, and eyes. Do not take internally.  
**Eye contact:** B4<sup>®</sup> Optimizer contains PEG which may be irritating to eyes. Before using this product wear protective glasses as well as covering the patient's eyes to protect from splashing material. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.  
**Skin contact:** B4<sup>®</sup> Optimizer contains PEG which can cause skin sensitization (allergic contact dermatitis) in susceptible individuals. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water after contact. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention.  
**Oral mucosa contact:** Avoid prolonged contact with oral soft tissues. If accidental contact occurs, flush mucosa with plenty of water and expectorate water. If sensitization of mucosa persists, seek medical attention immediately.
2. B4<sup>®</sup> Optimizer should not be used with patients who have a history of severe allergic reaction to any of the components.

### PRECAUTIONS

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
3. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
4. The B4<sup>®</sup> Optimizer bottle should be tightly closed immediately after use.
5. Use only in well ventilated areas.
6. Flammable: B4<sup>®</sup> Optimizer contains ethanol. Keep away from sources of ignition.

### STORAGE

B4<sup>®</sup> Optimizer should be kept out of direct sunlight and stored in a well ventilated place at room temperature not exceeding 25°C/77°F. Refrigerated storage is acceptable when not in use. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

### ADVERSE REACTIONS

1. Product may irritate the eyes and skin. Eye contact: irritation and possible corneal damage. Skin contact: irritation or possible allergic response. Reddish rashes may be seen on the skin. Mucous Membranes: inflammation, edema, sloughing. (See Warnings)
2. B4<sup>®</sup> Optimizer may cause pulpal effects. (See Contraindications)

### STEP-BY-STEP INSTRUCTIONS FOR USE

#### Pre-Impression Preparation Conditioner

1. **Preparation:** Complete tooth preparation, including core buildup or implant abutment/coping connection and/or retraction cord placement as necessary, in usual fashion.
2. **Cleaning:** Remove provisional restoration and residual temporary cement (if applicable). Clean preparation with a rubber cup and pumice or a cleaning paste such as NUPRO<sup>®</sup> Prophylaxis Paste. Wash thoroughly with water spray and air dry. Clean freshly instrumented preparation or implant component with water spray and then air dry.
3. **Cord placement:** If retraction cord is to be used, place per usual technique in accordance with manufacturer's instructions.
4. **Application of B4<sup>®</sup> Optimizer:**
  - 4.1 Dispense 2-3 drops of B4<sup>®</sup> Optimizer from the bottle into a clean dampen dish. Apply additional drops for multiple preparations, if needed. Replace cap promptly.
  - 4.2 Using a clean disposable brush, apply B4<sup>®</sup> Optimizer to the entire preparation/abutment/coping and sulcus, thoroughly (re)wetting the packed retraction cord, if applicable. Air-thin for a uniform film but do not rinse or dry applied B4<sup>®</sup> Optimizer.
  - 4.3 Remove cord, if applicable, and confirm hemostasis per usual technique. If additional chemical agents are used to establish hemostasis, thoroughly rinse preparation and sulcus and re-apply B4<sup>®</sup> Optimizer as outlined in section 4.2 once bleeding is controlled. Once hemostasis is confirmed, proceed immediately to impression as per Aquasil Ultra Impression Material's Directions for Use.

### CLEANING AND DISINFECTION

#### Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

To prevent B4<sup>®</sup> Optimizer bottles from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination. Repeated disinfection may damage label.

### LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY/MM"
2. The following numbers should be quoted in all correspondence:
  - Reorder Number
  - Lot number on bottle
  - Expiration date on bottle

Manufactured by:  
DENTSPLY Caulk  
38 West Clarke Avenue  
Milford, DE 19963 USA  
Tel.: 1-302-422-4511  
www.dentsply.com

DENTSPLY DeTrey GmbH  
De-Trey-Str. 1  
78467 Konstanz  
Germany  
Tel.: 49-7531-583-0  
www.dentsply.de

Distributed by:  
DENTSPLY Canada  
161 Vinyl Court  
Woodbridge, Ontario  
L4L 4A3 Canada

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